

DENTAL UNIT WATERLINE AND BOTTLE CLEANER

For the cleaning and control of microbial contamination in dental unit waterlines and bottles. For use in any dental unit.

DIRECTIONS FOR USE:

IT IS A VIOLATION OF FEDERAL LAW TO USE THIS PRODUCT IN A MANNER INCONSISTENT WITH ITS LABELING.

INITIAL CLEANING

***NOTE:** Each week allow Mint-A-Kleen to remain in the dental unit waterlines and bottle overnight and over the weekend.

Week 1: Use Mint-A-Kleen at the **end of every business day**.

Weeks 2-5: Use Mint-A-Kleen **twice per week** at the end of the business day.

After Week 5: Use Mint-A-Kleen **once per week**.

INITIAL CLEANING PURPOSE: Mint-A-Kleen preps the dental unit waterlines by gradually penetrating the microbial contamination during this 5 week period. This initial cleaning prevents copious amounts of contamination from clogging the waterlines.

ROUTINE WEEKLY CLEANING

- Use Mint-A-Kleen **once per week**.
- Use Mint-A-Kleen at the **end of the last business day of the week**.
- Leave Mint-A-Kleen in dental unit waterlines and bottle **overnight and over the weekend**.
- On Monday morning before dental procedures begin, Air Purge Waterlines again then Flush Waterlines with Water.


DIRECTIONS FOR WATERLINE AND BOTTLE CLEANING



STEP 1:
AIR PURGE WATERLINES



STEP 2:
ADD MINT-A-KLEEN*



STEP 3:
AIR PURGE WATERLINES



STEP 4:
FLUSH WATERLINES
WITH WATER

STEP 1: AIR PURGE: REMOVE DENTAL UNIT WATER FROM DENTAL UNIT WATERLINES AND BOTTLE

- After the last patient leaves at the end of the business day, turn dental unit **OFF**.
- Empty water bottle. Now attach this empty water bottle to dental unit.
- Turn dental unit **ON**.
- Air Purge waterlines by activating handpiece hoses and air/water syringes until water is eliminated.

STEP 2: ADD MINT-A-KLEEN

- Turn dental unit **OFF**. Remove water bottle; pour 4 ounces of Mint-A-Kleen into dental unit's empty water bottle.
- Place sterile 2x2 gauze on dental water bottle's opening. Turn dental water bottle upside down 3 times with gauze in place. This allows Mint-A-Kleen to coat this dental water bottle's interior walls.
- Using the same 2x2 gauze, wipe down the pick-up tube.
- Attach dental water bottle to dental unit.
- Turn dental unit **ON**. Activate handpiece hoses and air/water syringes until Mint-A-Kleen exits.
- Turn dental unit **OFF**. Allow Mint-A-Kleen to treat waterlines and bottle for at least 12 hours before proceeding to the next step.

STEP 3: (REPEAT STEP 1: AIR PURGE DENTAL UNIT WATERLINES)

- Turn dental unit **OFF**. Remove water bottle, dispose of any excess Mint-A-Kleen. Empty water bottle. Now attach this empty water bottle to dental unit.
- Turn dental unit **ON**.
- Air Purge Waterlines by activating handpiece hoses and air/water syringes until **NO** liquid comes out.

STEP 4: FLUSH DENTAL UNIT WATERLINES AND WATER BOTTLE

- Turn dental unit **OFF**. Remove water bottle from dental unit, dispose of any excess Mint-A-Kleen.
- Rinse dental water bottle with water. Fill water bottle with water source.
- Turn dental unit **ON**. Flush waterlines with water for 30 seconds to remove any excess Mint-A-Kleen.
- Begin Dental Procedures.

***NOTE:** For optimal cleaning, air purge and clean all instruments and tubing that contain water using the aforementioned protocol, including but not limited to: handpiece, air/water syringe, scaler, quick disconnect line.

Allow Mint-A-Kleen to remain in the dental waterlines while the dental office is closed for the weekend, holiday break or vacation. Dental Schools that use Mint-A-Kleen leave Mint-A-Kleen in their dental waterlines for up to 3 months during their summer breaks.

*****See Other Side of this Insert for Precautionary Statements, First Aid, Storage and Disposal.*****



PROVEN FORMULA FOR WATERLINE TREATMENT BY ANODIA SYSTEMS

See Other Side of this Insert for Directions for Use.



Mint-A-Kleen 1 PACK: MFR #: MAK-OZ-1 (16 oz Bottle)
Mint-A-Kleen 6 PACK: MFR #: MAK-OZ-6 (16 oz Bottles)
Mint-A-Kleen 10 PACK: MFR #: MAK-OZ (16 oz Bottles)

MANUFACTURED BY ANODIA SYSTEMS

EPA Reg. No. 85298-2

EPA Establishment No. 85298-KY-001

STORAGE AND DISPOSAL

Do not contaminate water, food or feed by storage and disposal.

Storage: store in a cool dark room away from heat or open flame.

Disposal: wastes resulting from the use of this product must be disposed of on site or at an approved waste disposal facility.

Container: non-refillable container. Do not reuse or refill container. Rinse with water and place in trash or offer for recycling if available.

PRECAUTIONARY STATEMENTS:

Hazards to Humans and Domestic Animals

CAUTION: Causes moderate eye irritation. Harmful if inhaled or swallowed. Avoid contact with eyes or clothing or breathing vapor. Remove and wash contaminated clothing before reuse. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum, using tobacco or using the toilet.

ACTIVE INGREDIENTS

Chlorhexidine Gluconate 0.12%

OTHER INGREDIENTS 99.88%

TOTAL 100.00%

CAUTION

KEEP OUT OF REACH OF CHILDREN

Net Contents: 16 FL OZ (473 ML)

FIRST AID

IF IN EYES: Hold eye open and rinse slowly and gently with water for 15-20 minutes. Remove contact lenses, if present, after the first 5 minutes, then continue rinsing. Call a poison control center or doctor for treatment advice.

IF INHALED: Move person to fresh air. If person is not breathing, call 911 or an ambulance, then give artificial respiration, preferably mouth-to-mouth if possible. Call a poison control center or doctor for further treatment advice.

IF SWALLOWED: Call a poison control center or doctor immediately for treatment advice. Have a person sip a glass of water if able to swallow. Do not induce vomiting unless told to by a poison control center or doctor. Do not give anything to an unconscious person. Have the product container or label with you when calling a poison control center or doctor or going for treatment.

PHYSICAL OR CHEMICAL HAZARDS

Combustible. Do not use or store near heat or open flame.

Mint-A-Kleen SDS available at www.Mint-A-Kleen.com

COPYRIGHT PROTECTED ©2009-2019 ANODIA SYSTEMS.

MINT-A-KLEEN Instructions and MINT-A-KLEEN Online Instructional Video are Registered with the United States Copyright Office. All Rights Reserved.

No part of this publication may be reproduced, stored in a retrieval system or transmitted in any form or by any means, mechanical, electronic, photocopying, recording, or otherwise without the prior written permission of the author, Anodia Systems.

MINT-A-KLEEN and MINTAKLEEN are Registered Trademarks of Anodia Systems. PATENTS PENDING.



BUILT FROM EXPERIENCE

109 Larrimore Lane | Danville, KY 40422 USA
P: 866.246.2548 | E: info@AnodiaSystems.com

www.AnodiaSystems.com