

# AirSTEALTH™

## Dental High-Speed Suction HVE

Patented by Anodia Systems, LLC

Item #: MAK-BLACK-HVE

### INDICATIONS FOR USE

AirSTEALTH™ is a Dental High-Speed Suction HVE featuring Anodia Systems' exclusive Hygienic & Leverless® Design Comprised of 3 Anodized, Autoclavable Parts: Base, Insert and Shroud. AirSTEALTH's Smooth Body Contour and Hidden Valve (Insert) Enables an ON / OFF Suction Feature by Sliding the Shroud Up and Down. For Cleaning & Sterilization Purposes, AirSTEALTH Assembles and Disassembles within Seconds!

### AirSTEALTH's 3 Parts:

#### AirSTEALTH includes:

1. Suction Cap
  - o Place the Enclosed Suction Cap Over the HVE's **Base** Opening when the **Shroud** and **Insert** are being Sterilized while the **Base** is still Connected to the Suction Hose. This Suction Cap Eliminates the Suction Noise.
2. One Assembled AirSTEALTH HVE Comprises of 3 Parts: Base, Shroud & Insert.



Assembled AirSTEALTH



Suction Cap  
Covers Base  
when not in use



Base



Shroud




Insert

## DIRECTIONS FOR USE

### Suction ON Position:

Turn Suction **ON**  
by Sliding  
**Shroud & Base Together.**



### Suction OFF Position:

Turn Suction **OFF**  
by Sliding  
**Shroud Away from Base.**

### ROUTINE CARE

1. After Each Patient, Disconnect AirSTEALTH's **Shroud** from AirSTEALTH's **Base**.
2. Leave AirSTEALTH's **Base** Connected to the Suction Hose Tubing. Place Suction Cap on **Base's** Opening.
3. Take AirSTEALTH's 2 Parts: **Shroud & Insert** to the Sterilization Area. Remove **Insert** from **Shroud**.
4. Install a Clean, Sterile AirSTEALTH **Insert**, **Shroud** into AirSTEALTH's **Base**, which is still Connected to the Suction Hose Tubing. Begin Dental Operative Procedures.

### IMPORTANT

- After Sterilization, Lubricate the **Insert's** 2 O-Rings and the **Base's** 2 O-Rings with a Lubricant such as Vaseline. Reassemble AirSTEALTH.
- Steam Autoclave, Dry Heat, Ethylene Oxide, and Gas or Chemical Vapor Process are All Safe to Use, as long as Temperatures Do Not Exceed 280 Degrees Fahrenheit (138 Celsius). **NOTE:** If you use the Chemical Vapor Process, AirSTEALTH must First be thoroughly Rinsed of All Cleaning Agents. Cleaning Agents left on AirSTEALTH may React with the Chemical Vapor Solutions to Form a Crystalline Residue, which can Clog AirSTEALTH or cause it to Stick.

**See other side for Assembly & Disassembly**

**ANODIA**®  
SYSTEMS  
*Dental Manufacturer*

[www.AnodiaSystems.com](http://www.AnodiaSystems.com)

See other side for more Product Info

## ASSEMBLY

**1**

Firmly Push Base (threaded end) Down into Suction Hose

**2**

Base Connected to Suction Hose

**3**

Slide Shroud Down onto Base

**4**

Position Insert's Hole with 2 Black O-Rings on Top

**5**

Drop Insert Down into Shroud

**6**

Firmly Push Insert Down to Base with Mirror Handle

Insert Clicks TWICE into Base

**7**

Suction ON

Slide Shroud Down to Base

**8**

Suction OFF

Slide Shroud Away from Base

## DISASSEMBLY

2 Steps:

**1**

Pull Shroud apart from Base

Base Stays Connected to Suction Hose

**2**

Push Insert through Shroud's Top Opening