

AirSTEALTH™

Dental Saliva Ejector

Patented by Anodia Systems, LLC

Item #: MAK-BLACK-SE

INDICATIONS FOR USE

AirSTEALTH™ is a Dental Saliva Ejector featuring Anodia Systems' exclusive Hygienic & Leverless® Design Comprised of 3 Anodized, Autoclavable Parts: Base, Insert and Shroud. AirSTEALTH's Smooth Body Contour and Hidden Valve enables an ON / OFF Suction Feature by Sliding the Shroud Up and Down. For Cleaning & Sterilization Purposes, AirSTEALTH Assembles and Disassembles within Seconds!



Saliva Ejector Suction Tip
[Removable & Autoclavable]



Assembled AirSTEALTH

AirSTEALTH's 3 Parts:



Base Shroud Insert



Suction Cap
Covers Base
when not in use




AirSTEALTH includes:

1. Suction Cap
 - o Place the Enclosed Suction Cap Over the Saliva Ejector's Base Opening when the Shroud and Insert are being Sterilized while the Base is still Connected to Suction Hose.
 - o This Suction Cap Eliminates the Suction Noise.
2. Saliva Ejector Suction Tip (Autoclavable)
 - o Place Saliva Ejector Tip on Top of Shroud During Dental Operative Procedures.
3. One Assembled AirSTEALTH Saliva Ejector Comprised of 3 Parts: Base, Shroud & Insert

DIRECTIONS FOR USE

Suction ON Position:

Turn Suction **ON**
by Sliding
Shroud Away from Base.



Suction OFF Position:

Turn Suction **OFF**
by Sliding
Shroud & Base together.

ROUTINE CARE

1. After Each Patient, Disconnect AirSTEALTH's **Shroud** from AirSTEALTH's **Base**.
2. Leave AirSTEALTH's **Base** Connected to the Suction Hose Tubing. Place Suction Cap on **Base's** Opening.
3. Take AirSTEALTH's 2 Parts: **Shroud & Insert** to the Sterilization Area. Remove **Insert & Tip** from **Shroud**.
4. Install a Clean, Sterile AirSTEALTH **Insert, Shroud** with **Tip** into AirSTEALTH's **Base**, which is still Connected to the Suction Hose Tubing. Begin Dental Operative Procedures.

IMPORTANT

- After Sterilization, Lubricate the **Insert's** 2 O-Rings and the **Base's** 2 O-Rings with a Lubricant such as Vaseline. Reassemble AirSTEALTH.
- Steam Autoclave, Dry Heat, Ethylene Oxide, and Gas or Chemical Vapor Process are All Safe to Use, as long as Temperatures Do Not Exceed 280 Degrees Fahrenheit (138 Celsius). **NOTE:** If you use the Chemical Vapor Process, AirSTEALTH must First be thoroughly Rinsed of All Cleaning Agents. Cleaning Agents left on AirSTEALTH may React with the Chemical Vapor Solutions to Form a Crystalline Residue, which can Clog AirSTEALTH or cause it to Stick.

See other side for Assembly & Disassembly

ANODIA
SYSTEMS
Dental Manufacturer

www.AnodiaSystems.com

AirSTEALTH™

See other side for more Product Info

ASSEMBLY

1 *Firmly Push*
Base (threaded end)
Down into
Suction Hose

2 Base
Connected to
Suction Hose

3 *Slide*
Shroud
Down onto
Base

4 *Firmly Push*
Insert
Down into
Shroud

5 *Push*
Insert
Down to
Base with
Mirror Handle

6 *Push*
Suction Tip
on
Shroud

7 Suction
ON

8 Suction
OFF

Slide
Shroud
Away from
Base

Slide
Shroud
Down to
Base

DISASSEMBLY

2 Steps:

1 *Pull*
Shroud
apart from
Base

2 *Push*
Insert
through
Shroud's
Top Opening

Base
Stays
Connected to
Suction Hose

ANODIA[®]
SYSTEMS
Dental Manufacturer

www.AnodiaSystems.com